PARSONS
ADL

Comfort Wash

Inflatable hair wash basin

CE

Necessary equipment:

• Comfort Wash • Towels • Water collection bucket • Shampoo • Water supply or pitchers of warm water

Instructions for use:

1. Inflate with mouth or pump. 2. Place Comfort Wash under patient's head. 3. Fold one towel and place between upper shoulders and Comfort Wash for comfort. 4. Position drain hose at least one foot higher than anticipated water level. *Note: A hose shut-off valve may be used as an added convenience for frequent use.* 5. Wet hair and shampoo as usual. 6. Rinse hair well and drain water into bucket. 7. Towel dry hair. 8. Wipe out Comfort Wash leaving inflated for several hours before deflating, folding and storing.

Limited One Year Warranty : Parsons ADL Inc. warrants to the original owner of this product that it is free from defects in material and workmanship for one year from date of purchase. We will repair or replace any part of the product, which in our opinion is defective, provided the product has not been abused, misused, altered or damaged after purchase. (This includes damage due to use of tools or harsh chemicals).

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