
Note: Prescribed treatments need to be very specific and individualized. Therefore, these instructions for using the Reciprocating Pulley Device are of a general nature only. For an exercise program specific to your needs, please consult your physician or physiotherapist.

Instructions:

- Place hook over the top of a door that can be closed tightly. For safety, place a warning card on the opposite doorknob so someone won't accidentally open the door while you are exercising.
- Adjust chair distance from door, or stand a distance from door so that the pulley is slightly behind your head and the length of rope is comfortable for your arm movement.
- As with any new exercise program, make sure you consult your physician or physiotherapist to ensure proper form when using the device.
- Complete exercises in a slow, smooth motion, being careful not to strain your arm muscles.
- If you have only one arm affected, you can use your other arm to assist in elevating the affected arm. Pull until you feel a small amount of tension and stretch in the affected arm. Hold this position for 4 to 7 seconds and then release. Relax the affected arm before the next repetition.
- As the affected arm gets stronger, you can apply a slight resistance with your other arm to increase the amount of muscular work done by the affected arm.
- Do only the assigned number of repetitions as prescribed by your medical practitioner.
- Performing your prescribed exercise program diligently (a typical program consists of exercising 3 times per week for 4 to 6 weeks) will help increase or maintain the active range of motion of your shoulders and arms, encourage normal joint mechanics and movement, assist in joint lubrication, and increase muscular strength. So become and stay dedicated to your program!

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