



SOLID SEAT INSERT

Designed to prevent wheelchair cushions from sagging when used in a sling seat. These inserts are made of 1/4" birch plywood for strength. All edges and corners are radiused to prevent damage to cushion covers. These inserts will help to promote better posture and therefore are more comfortable.

Support rigide pour coussin de fauteuil roulant

Conçu pour empêcher l'affaissement d'un coussin sur l'assise non rigide d'un fauteuil roulant. Fabriqué de contre-plaqué de bouleau, ce produit solide mesure 0,25 po (6,4 mm) d'épais. Les bords de ce support sont arrondis afin de ne pas endommager les housses de coussin. Surtout pratique pour promouvoir une meilleure posture et ajouter au confort.



Different sizes available - Plusieurs grandeurs offertes :

18C306 - 14 x 16 in/po (35.6 x 40.6 cm)	18C307 - 14 x 18 in/po (35.6 x 45.7 cm)
18C300 - 16 x 16 in/po (40.6 x 40.6 cm)	18C301 - 16 x 18 in/po (40.6 x 45.7 cm)
18C303 - 16 x 20 in/po (40.6 x 50.8 cm)	18C302 - 18 x 18 in/po (45.7 x 45.7 cm)
18C304 - 18 x 20 in/po (45.7 x 50.8 cm)	18C305 - 20 x 20 in/po (50.8 x 50.8 cm)
18C308 - 20 x 22 in/po (50.8 x 55.9 cm)	

Limited One Year Warranty

Parsons ADL warrants to the original owner of this product that it is free from defects in material and workmanship for one year from date of purchase. This warranty is offered provided the product has not been abused, misused, altered or damaged after purchase. (This includes damage due to the use of harsh chemicals.)

Garantie limitée d'un an

Parsons ADL garantit au propriétaire d'origine du présent produit que celui-ci est exempt de défauts de matériaux et de fabrication pour une durée d'une année à compter de la date d'achat. Cette garantie est offerte à condition que le produit n'ait pas été utilisé abusivement, mal utilisé, modifié ou endommagé après l'achat. (Ceci comprend les dommages résultant de l'utilisation de produits chimiques puissants.)

Parsons ADL Inc. 1986 Sideroad 15, Tottenham, Ontario L0G 1W0 Canada

Tel: 905.936.3580 - 1.800.263-1281 / Fax: 905.936.3585 - 1.888.225.5535

Email: custserv@parsonsadl.com

Aids to Daily Living - Accessoires à la vie quotidienne